# DIGITAL ELECTRONIC FEVER THERMOMETER

## INSTRUCTIONS FOR USE – Please read carefully

### Features:

- Easy to read digital liquid crystal display (LCD)
- Temperatures can be taken under arm (axillary use) or under the tongue (oral use)
- Auto shut-off after approximately 10 minutes of inactivity
- Waterproof for easy cleaning
- ABS resin body presents fewer hazards than mercury glass thermometers
- Low batter indicator
- Fast and accurate – 10 second read time after thermometer reaches a stable temperature reading

**Thermometer's Measurement Range:** 32.0~42.9°C or 90.0~109.9°F

**Resolution:** 0.1°C / °F

**Accuracy:** ±0.1°C, 35.5~42.0°C (±0.2°F, 95.9~107.6°F), ±0.2 C below 35.5°C or over 43.0°C (±0.4°F below 95.9°F or over 107.6°F)

**Battery:** LR 41-Alkaline button cell battery (0.0% Hg)

**Beeper:** Approximately 10 seconds once a stable temperature reading is reached

**Memory:** Built-in memory enabling display of previously measured value

**°F/°C Change:** While unit is off, press & hold ON/Off button for approximately 4 seconds to toggle between °F/°C. Release the button once desired setting is established

**Working Conditions:** Relative Humidity: 15% ~ 95%, Temperature: 5~40°C

**Storage Conditions:** Relative Humidity: 15% ~ 95%, Temperature: -20~55°C

### Battery Replacement:


1. Replace battery when "⚡" appears in lower right corner of LCD display
2. Carefully pull battery cover off (located at opposite end of probe)
3. Gently pull out plastic circuit board with battery chamber approximately 1 cm. Pulling the battery circuit board out too fast or more than 1cm will damage the internal wiring
4. Use a small instrument like a toothpick to remove battery. Avoid using a sharp metal object
5. Place a battery in the chamber with positive side (+) facing up
6. Slide battery chamber back into place and put the cover on
7. Warning: discard only battery carefully. Do not dispose of the battery in a fire. It may explode

### Caution:

- Read instructions completely before using digital thermometer
- Keep out of reach of children and pets
- Thermometer cap and battery may be fatal if swallowed. If this occurs, seek medical attention immediately
- Do not use thermometer in ear. Designed for oral and axillary readings only
- Do not place thermometer battery near extreme heat
- Note: use of probe cover may result in 0.1°C (0.2°F) discrepancy from actual temperature
- Normal body temperature varies with age, time of day, physical activity, hormonal and environmental factors
- Temperature in the late afternoon or evening will usually be 0.6°C (1.0°F) higher than the morning temperature
- It is recommended to check the performance of the thermometer every 2 years. Do not disassemble
- Performance of thermometer may be affected if the following occurs:
  - Storage or operation outside stated temperature and humidity range
  - Mechanical shock
  - Electromagnetic interferences (ie. microwave, cell phones)

### Instructions for use:

1. Disinfect the probe with isopropyl (rubbing) alcohol. **DO NOT** boil probe.
2. Press the on/off button next to LCD display. A tone will sound as the screen displays **|||||** followed by the last recorded temperature.
3. The display will flash **"Lo"** (°F or °C). The last recorded temperature is displayed before the screen display goes back to **"Lo"** (°F or °C). Note: if the measured or room temperature is less than 90°F or 32°C, the LCD will display **Lo**, if the measured or room temperature is higher than 109.9°F or 42.9°C, the LCD will display **Hi**. Lo or Hi does not indicate defective battery/thermometer.

4. Position thermometer in desired location (mouth or armpit).
    - Oral use:** Place thermometer under the tongue as indicated with a check mark. Close mouth and breathe through the nose. Normal temperature should be 36.3°C to 37.2°C (97.3°F to 99.1°F).
    - Armpit use:** Wipe armpit dry. Please probe in armpit and keep arm pressed firmly at side. Normal temperature should be 35.8°C to 36.8°C (96.4°F to 98.2°F).
  5. The degree sign continues flashing on the LCD while the thermometer is reading your temperature. 10 seconds after thermometer reaches a stable temperature reading, the flashing will stop, and an alarm will beep. The measured reading will appear on the LCD. The displayed reading will not change when the thermometer is removed from its testing position.
    - Normal alarm:** Steady beeps Bi---Bi---Bi.
    - High temperature alarms:** Rapid and lengthy beeps Bi-Bi-Bi..... Bi-Bi-Bi..... Bi-Bi-Bi.
  6. To prolong battery life, turn unit off after testing. If no action is taken, unit will automatically shut off after approximately 10 minutes.
  7. Disinfect the thermometer and store in the case provided.
- For optimal product performance, it is recommended to replace thermometers after 5 years of usage.

**Questions? 1.888.475.2500 | [www.pharmasystems.com](http://www.pharmasystems.com)**

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Manufactured for / Produit pour: PharmaSystems Inc.