

## Conical Glass Graduate, 5 mL



Item #: 3360  
UPC: 063636 833602

- 5 mL conical glass cylinder
- Marked in increments of 1 mL
- Bell shaped glass graduated cylinders, with easy-pour spouts
- Bell shape has a controlled diameter – recommend for accuracy
- White enamel graduations, single scale
- Volumes calibrated to deliver (TD)
- Can be autoclaved

\* Discount pricing available. Log in to see pricing and details.

### Available Sizes:

Item	Volume	Shape	Material	Increments
3360	5 mL	Bell	Borosilicate	1 mL
3361	10 mL	Bell	Borosilicate	0.5 mL
3362	25 mL	Bell	Borosilicate	1 mL
3019	50 mL	Conical	Flint	5 mL
3021	100 mL	Conical	Flint	5-25 mL
3022	250 mL	Conical	Flint	10-50 mL
3023	500 mL	Conical	Flint	10 mL
3363	1000 mL	Bell	Borosilicate	20 mL

